


Eagle 620 / 625

Patient lift



The patient lifts Eagle 620 and Eagle 625 combines functionality, safety, and reliability. The qualitative wheels combined with a narrow total width allow for easy transfers, even in narrow spaces.



VERMEIREN
we care for you

Slings

Different lifting slings are available for our patient lifts. Whether you opt for the standard, hygienic, or comfort version with upper body support, our lifting straps offer a solution suitable for everyone's needs for support. The different sizes are indicated with color on the slings.

- Medium (M) yellow
- Large (L) blue
- Extra Large (XL) black

Patient weight:

M (yellow) 40 - 90 kg
L (blue) 70 - 120 kg
XL (black) 100 - 200 kg

Colours

C29
grey



Technical data Eagle 620

Total length	1115 mm
Total width (closed) / (opened)	550 / 1035 mm
Max. total height	1740 mm
Weight	37 kg
Maximum user weight	150 kg
Additional information	CE

Products designed to suit you


VERMEIREN always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.

SLING HYGIENE

Available in

M - L - XL - standard or netted structure

	A	B	C	D
M	780	250	600	300
L	840	260	700	360
XL	910	270	820	440




SLING COMFORT

Available in

M - L - XL - standard or netted structure

	A	B	C	D
M	780	400	450	300
L	840	450	500	360
XL	910	510	560	440




SLING STANDARD

Available in

M - L - XL - standard or netted structure

	A	B	C	D
M	780	400	450	300
L	840	450	500	360
X	910	510	560	440



Eagle 625

Total length	1290 mm
Total width (closed) / (opened)	680 / 1225 mm
Max. total height	1840 mm
Weight	42 kg
Maximum user weight	175 kg
Additional information	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road
HR4 9NS Hereford - United Kingdom
Company registration number 16075613

www.vermeiren.co.uk
info@vermeiren.co.uk



For more information, contact your local dealer

VERMEIREN
we care for you