



Albatros

Active lift

The Albatros is a compact patient hoist with easy usage.
Foldable in a few seconds for easy storage and transport.



VERMEIREN
we care for you



SLING SIZES	
Available in M - XL (standard or netted structure)	
M	XL
40-120 KG	100-150 KG

Products designed to suit you

VERMEIREN always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Albatros** complies with the highest norms, laid down by certified test institutes in most European countries.

The ACTIVE sling is available in 2 sizes and is washable at 60° Celsius. The medium size supports patients from 40 up to 120 kg, whilst the XL size supports those from 100 up to 150 kg.

Colours

C29
grey



Technical data

■ Total width (closed)	584 mm
■ Total width (opened)	945 mm
■ Total height min.	1040 mm
■ Total height max.	1739 mm
■ Total length	1046 mm
■ Folded height	500 mm



■	
■	
■ Leg height / Frame height	117 mm
■ Weight	44,85 kg
■ Maximum user weight	150 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road
HR4 9NS Hereford - United Kingdom
Company registration number 16075613

www.vermeiren.co.uk
info@vermeiren.co.uk



For more information, contact your local dealer

VERMEIREN
we care for you