

Bobby Evo

Compact, lightweight and practical transfer wheelchair

The Bobby Evo is your ideal compact, foldable transferchair. Light and easy to fold in seconds.

- Swing away armrests for easy lateral transfers
- Reduced total width for optimal passage in tight spaces (+7cm compared to seat width)



VERMEIREN
we care for you

Features

- Padded seat and back upholstery
- Storage bag on the back of the backrest
- Padded swing away armrests (B03)
- Semi-foldable backrest
- Ergonomic push handles
- Brakes for the accompanying person
- 2 manual brakes with long handle
- Angle adjustable footplates
- Adjustable ankle straps
- Safety belt

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalised so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Bobby Evo** complies with the highest norms, laid down by certified test institutes in most European countries.



Technical data

Total width	460 - 490 - 550 mm
Total height	900 mm
Total length	1010 mm
Folded width	260 mm
Seat width	390 - 420 - 480 mm
Height armrests	200 mm

Seat height	490 mm
Height backrest	440 mm
Seat depth	420 mm
Weight	12,25 * kg
Maximum user weight	115 kg
Additional info	CE

* configuration with B03 armrest
Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible
Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road
HR4 9NS Hereford - United Kingdom
Company registration number 16075613

www.vermeiren.co.uk
info@vermeiren.co.uk



For more information, contact your local dealer