

Jazz S50

Standard wheelchair

The wheelchair **Jazz S50** is an excellent combination between technicality, comfort and **VERMEIREN** quality. Often seen as a standard wheelchair but with many adjustability options.

Ergonomic handgrips, easy to clean nylon upholstery, adjustable arm- and legrests and detachable rear wheels.



VERMEIREN
we care for you

Options

The **Jazz S50** can be equipped as transit model (T30) or extended with several options such as travelwheels, anti-tipping or a B05 armrest. Please consult your dealer or the orderform for a complete and actual list.



B78 - Anti-tipping wheels



B86 - Travelwheels



B05 - Adjustable armrest



B69 - Foldable backrest

Colour(s)

C86
carbon
grey



B75 - Brakes

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Jazz S50** complies with the highest norms, laid down by certified test institutes in most European countries.



T30 - Transit version

Technical data

Total width	565 / 595 / 615 / 635 / 655 / 675 mm
Folded width	290 mm
Seat with	390 / 420 / 440 / 460 / 480 / 500 mm
Height armrests	220 mm
Seat height	500 mm
Height backrest	420 mm

Total height	935 mm
Seat depth	450 mm
Total length	1120 mm
Weight	16,6 kg
Users weight	130 kg
Additional info	 

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road
HR4 9NS Hereford - United Kingdom
Company registration number 16075613

www.vermeiren.co.uk
info@vermeiren.co.uk



For more information, contact your local dealer

VERMEIREN
we care for you