



B L U M I L



distributed by





# Patented Connecting System



The Blumil add-ons have one incredibly unique feature - a patented connecting system. There is no need to attach any additional parts to your own wheelchair. You just press a button and an the electric system does the rest for you. 5 seconds later you are ready to GO with absolutely no assistance from anyone. You can be fully independent, no matter how much strength you have in your hands.





# B L U M I L -GO



## Speed modes

20 km/h

## Range

Up to 30 km

## Power

300 W

## Removable Battery

Li-ion 280 Wh

## Total weight

13 kg

## Wheel diameter

8,5"

## Lights

Built in LED lights

## Tire

Pneumatic, tubeless,  
puncture proof

## Connecting mechanism

Fully electric

## Certification

CE, medical device in EU

## Climbing capability

Up to 10 %

## Digital LCD-display

Current speed and speed  
mode, charge level, total  
distance

## Braking system

Drum brake and  
electromagnetic brake





# B L U M I L *SPORT*



## Speed modes

25 km/h + reverse

## Range

Up to 60 km

## Power

500 W

## Removable Battery

Li-ion 14,5 Ah

## Total weight

18 kg

## Wheel diameter

12"

## Lights

Powerful LEDx4

## Tire

Tube tire

## Certification

CE, medical device in EU

## Climbing capability

Up to 14 %

## Digital LCD-display

Current speed and speed mode, charge level, total distance, light indicator, average speed, top speed

## Braking system

Disc brake





CERTIFICATION



### CERTIFICATION

Blumil GO is a certified product. As a medical device it meets all the requirements stipulated by EU law. The CE mark itself is a determinant of certain safety and quality standards and provides information and proof to all consumers that we live up to strict requirements.



MADE IN EUROPE

### MADE IN THE EU

Blumil GO is of very high quality. It is manufactured in the EU to the highest standards.

distributed by

### VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre  
Faraday Road  
HR4 9NS Hereford  
United Kingdom  
Company registration number 16075613

[www.vermeiren.co.uk](http://www.vermeiren.co.uk)  
[info@vermeiren.co.uk](mailto:info@vermeiren.co.uk)



Contact your local supplier for more information.

**VERMEIREN**  
we care for you