



708D

Standard manual wheelchair

Aesthetics, reliability and functionality. All removable parts such as footrests, armrests, wheels and cushions can be replaced thanks to the modular system. The wide range of accessories and options allows the wheelchair to be optimally adapted to the user. Well known and recognized for its sturdiness and quality, the 708D makes an ideal asset for a rental fleet or as a durable choice to own for a very long time.



VERMEIREN
we care for you

Features



Aluminium footplates complete with rubber protection wheels



Easy to fold, easy to transport



Hemiplegic (HEM2) version. Double rims to propel and to steer



Amputee (AMP) version. Position of the rear wheel has shifted to the back

- 3 types of folding and removable armrests with long and comfortable armrests
- Wide range of options and accessories
- Armrest locking system is easy to access and to operate
- Various seat and backrest options
- Available in 6 seat widths
- 4 types of footrests and 5 types of leg rests

Colours



C29
grey

Technical data

■ Total width	560 / 590 / 610 / 630 / 650 / 670 mm
■ Total height	910 mm
■ Total length	1109 mm
■ Folded width	280 mm
■ Seat width	390 / 420 / 440 / 460 / 480 / 500 mm
■ Height armrests	225 mm

■ Seat height	490 mm
■ Height backrest	420 mm
■ Seat depth	430 mm
■ Weight	18 kg
■ Maximum user weight	120 kg
■ Additional info	

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **708D** complies with the highest norms, laid down by certified test institutes in most European countries.

VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road
HR4 9NS Hereford - United Kingdom
Company registration number 16075613

www.vermeiren.co.uk
info@vermeiren.co.uk



For more information, contact your local dealer

VERMEIREN
we care for you