



V500 Light

**OUR LIGHTEST
STANDARD WHEELCHAIR**



VERMEIREN
we care for you



Our lightest

Lighter than the V500, the V500 Light is our lightest standard wheelchair.

Whether for temporary or permanent use, it improves the user's daily life and makes transport and transfers easier.

Highly configurable

As part of the V-series, the V500 Light offers a wide range of accessories and options to create and modify your V500 Light.



Adjustable seat height (4 positions)

Adjustable seat depth: 42 to 48 cm



**Toolfree removable
armrests**



**In and outward folding
footrests
Removable footrests**



**Folding brakes levers
for easy transfers**



**Foldable for easy
transport or storage**



**Height adjustable push handles
and backrest**



Adjustable seat angle: 0° - 10°



**Adjustable armrests
in height and in depth**



V500 Light

Available colours



The **V500 Light** offers a wide range of accessories and options. Consult our online orderform to get a complete overview.

Technical data

Total width	590 - 620 - 640 - 660 - 680 - 700 mm
Total height	860 - 920 mm
Total length	975 - 1035 mm
Folded width	310 mm
Seat width	390 - 420 - 440 - 460 - 480 - 500 mm
Height armrests	150 - 210 mm

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalised so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Seat height	470 - 530 mm
Height backrest	400 - 420 - 440 - 460 - 480 mm
Seat depth	420 - 480 mm
Weight (without cushion)	15,5 kg
Maximum user weight	130 kg
Additional info	

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road
HR4 9NS Hereford - United Kingdom
Company registration number 16075613

www.vermeiren.co.uk
info@vermeiren.co.uk

For more information, contact your local dealer

