



D200 SPLIT

Lightweight manual wheelchair
with half folding backrest

The **D200 Split** is a lightweight aluminum manual wheelchair featuring quick release wheels and a comfort cushion and back support overlay in breathable spacer fabric. Attendant brakes come as standard along with a half folding back for ease of transportation.



VERMEIREN
we care for you

Custom-Designed Products

Our core principle is simple:

Your wheelchair should adapt to you, not the other way around. That's why our wheelchairs are evolutionary, offering the right response to your specific needs.

All our wheelchairs are adjustable, come with numerous options, and can be fully customised, ensuring there is a perfect wheelchair for every person. This is the result of our continuous research and experience since 1957.

Guaranteed Quality

In terms of quality, the **D200 Split** meets the highest standards, established by renowned certified testing institutes worldwide. For added peace of mind the frame comes with a 5 year warranty against manufacturing defects.

Colour

C86

Carbon Grey



Half folding back for easier transportation



Removable comfort cushion in breathable material



Attendant operated brakes with locking function

Options



B31 Crutch holder



B20 Safety belt



B78 Anti-tippers (Pair)



B66MA/B66OMA
A wide range of hemiplegic arm pads

Technical data

Overall width	580/600/620/640/660/680/700/720 mm
Overall height	880 - 940 mm
Overall length	1010 mm
Folded width	300 mm
Seat width	380/400/420/440/460/480/500/520 mm
Armrest height	220-240 mm

Seat height	470 - 530 mm
Backrest height	420 mm
Seat depth	430 mm
Weight	16 kg
Maximum user weight	130 kg
Additional information	 

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road
HR4 9NS Hereford - United Kingdom
Company registration number 16075613

www.vermeiren.co.uk
info@vermeiren.co.uk



For more information, contact your local dealer

VERMEIREN
we care for you