

Interval XXL

Care bed - 120 / 140 cm

The 3-functions Interval XXL bed has electrical adjustments and is available in 120 and 140 cm. Delivered in a white or beige frame with an extra strong, grid mattress base. The legpart comes with an electrically adjustable knee-buckle whereas the footpart can be adjusted manually. The complete bed can be ordered with wooden or metal side rails.



VERMEIREN
we care for you

Comfortable and practical

The Interval XXL bed is equipped with an electrically adjustable head part and knee buckle. With the simple touch of a remote button, you can position the bed to read, watching TV or dine. The leg rest in its turn can also be adjusted to the most comfortable position.

This variable height bed is equipped with 4 wheels with separate brakes, ensuring its stability at all times.

Aesthetic and easy to maintain

The foot and head board are an integral part of this XXL bed. The warm colour of the wood makes this bed an elegant piece of furniture that will fit in almost any interior.

Frame colors



*Easy-to-use remote
with safety locks*



*Metal mattress holders
on all sides*



*Manual adjustment
of the leg part*



*Easy handling knobs
to lower or raise
the side rails*

Technical data (120 cm mattress bed)

Total width / length / height	1330 x 2080 x 860 - 1295 mm
Total height of mattress support	370 - 805 mm
Mattress spec.	1200 x 2000 x 120 mm
Weight	124 kg
Max. load (incl. patient, mattress and options)	275 kg
Additional info	CE

Technical data (140 cm mattress bed)

Total width / length / height	1530 x 2080 x 860 - 1295 mm
Total height of mattress support	370 - 805 mm
Mattress spec.	1400 x 2000 x 140 mm
Weight	139 kg
Max. load (incl. patient, mattress and options)	275 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road
HR4 9NS Hereford - United Kingdom
Company registration number 16075613

www.vermeiren.co.uk
info@vermeiren.co.uk



For more information, contact your local dealer

VERMEIREN
we care for you