



V500 ACTIVE

Our heritage of the V-series
made active



VERMEIREN
we care for you



LIGHTNESS: FROM 11,9 KG

Through our efforts, we have successfully crafted an exceptionally lightweight model, weighing in at an impressive 11.9 kg. This remarkable outcome showcases our commitment to innovation and delivering a wheelchair that not only looks dynamic but also offers enhanced mobility and ease of use.

V500 ACTIVE


Introducing the **V500 Active wheelchair**, built upon the successful concept of the V500 family. This cutting-edge model brings a plethora of options and accessories, catering to diverse needs, including the most advanced requirements.

Designed for ultimate user personalisation, the V500 Active offers high configurability, ensuring a perfect fit for every individual. Embracing an active and modern aesthetic, it appeals to those who seek a stylish and dynamic appearance.


Experience the ideal fusion of comfort and adaptability found in a standard wheelchair, seamlessly combined with the sleekness and agility of an active folding chair. The V500 Active strikes the perfect balance, offering an unparalleled mobility solution that empowers users in their daily activities.



**ADJUSTABLE
SEAT ANGLE:**
0° - 10°



**ADJUSTABLE
SEAT HEIGHT**
44* - 47 - 50 - 53 CM



**ADJUSTABLE
SEAT DEPTH:**
42 TO 48 CM



**REMOVABLE,
IN AND OUTWARD
FOLDING FOOTRESTS**

**requires 22" rear wheels*



BACKREST HEIGHT

- 30 cm
- 35 cm
- 40 cm

PUSH HANDLES 5 SOLUTIONS

- Push handles - fixed (12 cm)
- Push handles - fixed short (8 cm)
- Foldable push handles
- Height adjustable fixed handles, with clamps
- Height adjustable foldable handles, with clamps



ARMRESTS

Armrests come in a spectrum of options, ranging from sleek aluminium side guards to convenient, adjustable and removable armrests designed for utmost comfort

- B82S - Clothing protectors (standard)
- B82 - Mud guards on clothing protector
- Bi5 - Removable height adjustable armrest
- Bi5A - Removable height adjustable curved mudguard
- Bi5SA - Removable height adjustable straight mudguard

FOOTRESTS

- Footrests and leg rests are offered in a diverse array of solutions, catering to various user requirements and adapting to their changing needs over time.

BRAKES

Push

Sport



WHEELS

- 11 solutions of rear wheels
- 21 solutions of front wheels
- To find your right balance between weight, comfort and budget.





AVAILABLE COLOURS

Embrace a spectrum of colours! Each shade we've carefully curated not only reflects our mood but also defines our way of life.

Unlock the opportunity to choose from a variety of hues, allowing you to personalise your wheelchair and truly embody your uniqueness.



C86 Carbon grey



C79 Light blue



C89 Copper



C85 Bright green



C76 White

Technical data

Total width	590 - 620 - 640 - 660 - 680 - 700 mm
Total height when folded (depends on seat height)	740 - 840 mm
Total length	975 - 1035 mm
Folded width	310 mm
Seat width	390 - 420 - 440 - 460 - 480 - 500 mm
Heaviest part	7,5 kg

Seat height	440 requires 22" rear wheels - 470 - 500 - 530 mm
Height backrest	300 - 350 - 400 mm
Seat depth	420 - 440 - 460 - 480 mm
Weight (without cushion)	From 11,9 kg
Maximum user weight	130 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road
HR4 9NS Hereford - United Kingdom
Company registration number 16075613

www.vermeiren.co.uk
info@vermeiren.co.uk



For more information, contact your local dealer

VERMEIREN
we care for you