



# D200

## Lightweight and functional

The D200 is a **lightweight manual wheelchair made of aluminium**. Foldable for easy transport and storage without comprimizing on adjustment possibilities. The new design and carbon colour give it a modern and agreeable look for all.



**VERMEIREN**  
we care for you

## Products designed to suit you


Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

## Quality guaranteed

Regarding quality, the **D200** complies with the highest norms, laid down by certified test institutes in most European countries.

## Colours



Rear wheels with quick-release system



Removable swing-away leg rests (in and outwards)



Option B78 - Anti-tipping wheels



Option BZ8 - articulated legrest self-correcting in length



Height adjustable push handles



Armrests adjustable in height and depth to suit specific need



Adjustable front wheels



Two tipping pieces

## Technical data

Total width	580/600/620/640/660/680/700/720 mm
Total length	1010 mm
Total height	880 - 940 mm
Folded width	300 mm
Seat width	380/400/420/440/460/480/500/520 mm
Height armrests	220 - 240 mm

Seat height	470 - 530 mm
Height backrest	420 mm
Seat depth	430 mm
Weight	16 kg
Maximum user weight	130 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

## VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road  
HR4 9NS Hereford - United Kingdom  
Company registration number 16075613

[www.vermeiren.co.uk](http://www.vermeiren.co.uk)  
[info@vermeiren.co.uk](mailto:info@vermeiren.co.uk)



For more information, contact your local dealer

**VERMEIREN**  
we care for you