



one

Compact indoor scooter
Perfect for daily trips
Adapted to your style

Introducing **One**, the fusion of compactness and style in a scooter designed for both indoor and adapted outdoor use. Its lightweight and compact design make it the ultimate companion for your daily journeys, effortlessly navigating through any urban setting. Embrace the freedom and convenience that **One** brings to your life, elevating your mobility experience to new heights. Experience the joy of owning a scooter that perfectly complements your style and adapts to your dynamic lifestyle.

one





Experience the power of **One = 4**

Adapt your scooter to your mood and your outfit

Scooter **One** comes with 4 interchangeable colours included. Express yourself effortlessly by easily swapping out the colours to create the perfect look for any occasion, in a few seconds. One scooter, four amazing styles - it's the ultimate equation for individuality and personalisation.

No tool needed





Easy to use

Designed with a 360° pivoting seat and foldaway armrest to ensures effortless accessibility, allowing users to easily get in and out of the scooter.



User-friendly controller with a battery status LED indicator, making it convenient and straightforward to operate.





Adjustable

Adjustable steering column, enabling users to effortlessly find their optimal driving position for maximum comfort.



Compact and efficient

Fully dismantlable in a matter of seconds the **One** scooter allows for effortless storage and travel, making it exceptionally convenient for users.



Fits in a small car

With an impressive minimum turning diameter of 115 cm, the **One** scooter excels as a compact indoor mobility solution, ideal for daily trips to the supermarket, restaurant visits, or manoeuvring through small spaces with ease.





one

Scooter

4 VERSIONS



one

Your **everyday** companion

8
km/h

16 km
driving range up to

15 Ah battery (24 V) | Solid PU tires



one



one air+

Go further with **comfort**

8
km/h



30 km



driving range up to

25 Ah battery (24 V) | Air tires



one air+

one^{LED} air+

See and be seen

8
km/h



driving range up to

30 km



25 Ah battery (24 V) | Lighting system



one^{LED} air+



one ^{SPRING} air+

Smoother rides, everytime

8
km/h

26 km
driving range up to

20 Ah battery (24 V) | Shock absorbers



one ^{SPRING} air+



	Lighting system	Shock absorbers	Tires	Speed	Battery type
one	X	X	PU	Up to 8 km/h	Gel
one 	X	X	Air	Up to 8 km/h	Gel
one ^{LED} 	✓	X	Air	Up to 8 km/h	Gel
one ^{SPRING} 	X	✓	Air	Up to 8 km/h	Lithium

	Battery	Driving range	Battery weight	Total weight
one	15 Ah (gel)	Up to 16 km	9 kg	44,9 kg
one 	25 Ah (gel)	Up to 30 km	14 kg	50 kg
one ^{LED} 	25 Ah (gel)	Up to 30 km	14 kg	50 kg
one ^{SPRING} 	20 AH (Lithium)	Up to 26 km	4,25 kg	42,35 kg



Technical information

■ Total width	490 mm
■ Total height	880 mm
■ Total lenght	1030 mm
■ Seat height	525-570 mm
■ Turning radius	1150 mm
■ Battery	24V - 15 Ah

■ Driving range	16 km*
■ Motors	270 W
■ Speed	8 km/h*
■ Weight	44,9 kg
■ Maximum user weight	136 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible Subject to technical alterations. Images for illustration purposes only. Certain features may be optional. * Theoretical range, calculated in accordance with EU standards under ideal conditions. Speeds and range may deviate (approx. 15%).



Technical information

■ Total width	490 mm
■ Total height	880 mm
■ Total lenght	1030 mm
■ Seat height	525-570 mm
■ Turning radius	1150 mm
■ Battery	24V - 25 Ah

■ Driving range	30 km*
■ Motors	270 W
■ Speed	8 km/h*
■ Weight	50 kg
■ Maximum user weight	136 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible Subject to technical alterations. Images for illustration purposes only. Certain features may be optional. * Theoretical range, calculated in accordance with EU standards under ideal conditions. Speeds and range may deviate (approx. 15%).



Technical information

■ Total width	490 mm
■ Total height	880 mm
■ Total lenght	1030 mm
■ Seat height	525-570 mm
■ Turning radius	1150 mm
■ Battery	24V - 25 Ah

■ Driving range	30 km*
■ Motors	270 W
■ Speed	8 km/h*
■ Weight	50 kg
■ Maximum user weight	136 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible Subject to technical alterations. Images for illustration purposes only. Certain features may be optional. * Theoretical range, calculated in accordance with EU standards under ideal conditions. Speeds and range may deviate (approx. 15%).



Technical information

■ Total width	490 mm
■ Total height	880 mm
■ Total lenght	1030 mm
■ Seat height	525-570 mm
■ Turning radius	1150 mm
■ Battery	24V - 20 Ah

■ Driving range	26 km*
■ Motors	270 W
■ Speed	8 km/h*
■ Weight	42,35 kg
■ Maximum user weight	136 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible Subject to technical alterations. Images for illustration purposes only. Certain features may be optional. * Theoretical range, calculated in accordance with EU standards under ideal conditions. Speeds and range may deviate (approx. 15%).

one



R.E.: **VERMEIREN GROUP**, Vermeirenplein 1/15 - 2920 Kalmthout - Belgium - BL One UK 2025 – 12v1.0

VERMEIREN
we care for you