

MAX. USER WEIGHT
250 KG

Forest 3+

bariatric wheelchair

The Forest 3+ is an electric wheelchair specially designed for large seat widths and for a higher user weight up to 250 kg. Complete with the latest electronic controller and a wide range of options.



VERMEIREN
we care for you

Adaptability

Compact, solid and efficient. The FOREST 3+ is the ideal companion for users who are looking for a compact wheelchair for indoors and at the same time a high-performance model for outdoors. Strength, comfort and adaptability are the main features of the FOREST 3+.

The front oscillating axle ensures that all four wheels contact the ground and ensures a smooth ride. The optional V-STAB technology allows the user to ride in a straight line, even on uneven terrain, without having to constantly compensate. The user will appreciate the comfort of the FOREST 3+ with its preformed padding and individual adjustment possibilities which make this seat highly adaptable.



Colours (matte and glossy colours)



We offer more than 32 different colours. For more information, please check our website.

Options



SE09 - Additional controller



SE04 Chin controller



SE55 - Reinforced joystick



L58 - Multi-adjustable headrest



SE83 - Mini joystick



B66MAL Hemiplegic armrest

Technical data

Total width	705 / 720 / 760 / 820 / 870 mm
Total length	1145 - 1163 mm
Total height	1100 mm
Seat height	450 mm
Seat width	500 / 550 / 600 / 650 / 700 mm
Seat depth	450 - 550 mm

Batteries	2 x 12V / 80Ah - 350 W
Speed*	6 km/h / 10 km/h
Driving range*	49 km
Weight	134,35 kg
Maximum user weight	250 kg
Additional information	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road
HR4 9NS Hereford - United Kingdom
Company registration number 16075613

www.vermeiren.co.uk
info@vermeiren.co.uk



For more information, contact your local dealer