



V100 XL

Maximised duties

The V100 XL is the entry-level model of the V-series, has a steel frame and is perfect for all-around use.



VERMEIREN
we care for you

Heavy duty wheelchair

The **V100 XL** is a robust and sturdy steel wheelchair with a user weight up to 170 kg with a wide range of wheelchair personalisation options. 100% customised to you and your needs. The double cross offers additional stability and 3 seat heights, 3 armrests and 3 upholstery options are creating an extremely versatile wheelchair.

Colours

C59

- ▶ 3 Seat widths
- ▶ 3 seat heights
- ▶ 3 armrest positions
- ▶ Maximum user weight: 170 kg
- ▶ Different options

Technical data

■ Total width	740 / 770 / 810 mm
■ Total height	890 - 950 mm
■ Total length	1060 mm
■ Folded width	375 mm
■ Seat width	530 / 560 / 600 mm
■ Height armrests	220 - 240 mm

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalised so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **V100 XL** complies with the highest norms, laid down by certified test institutes in most European countries.



■ Seat height	470 - 530 mm
■ Height backrest	380 mm
■ Seat depth	490 mm
■ Weight	21,4 kg
■ Maximum user weight	170 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road
HR4 9NS Hereford - United Kingdom
Company registration number 16075613

www.vermeiren.co.uk
info@vermeiren.co.uk



For more information, contact your local dealer

VERMEIREN
we care for you