

FALCON

Patient Lift



VERMEIREN
we care for you



Falcon, a safe and secure choice

A perfect mix between safety, security, and technology. Falcon allows for easy and comfortable transferring of patients. The two brakes assure a locked position for

ultimate safety when lifting or lowering. With its balanced design, maneuvering becomes smooth sailing for both patient and attendant.



Neoprene protection on spreader bar and 4 hooks for better stability



Ergonomic push handles



Controller with battery status indicator



Only 48 seconds from lowest to highest position



Battery status indicator



Audible warning when the battery is very low

Quality guaranteed

Regarding quality, **Falcon** complies with the highest norms, laid down by certified test institutes in most European countries.

Products designed to suit you

VERMEIREN always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.



Dismountable for easy transport and storage


Optional slings

REGULAR

Available in

M - L - XL (sizes in mm)

	A	B	C	D
M	780	400	450	300
L	840	450	500	360
XL	910	510	560	440




COMFORT

Available in

M - L - XL (sizes in mm)

	A	B	C	D	E
M	780	400	450	300	390
L	840	450	500	360	410
XL	910	510	560	440	420




HYGIENIC

Available in

M - L - XL (sizes in mm)


	A	B	C	D
M	780	250	600	300
L	840	260	700	360
XL	910	270	820	440



Optional slings - fabric



Polyethylene



Net Polyester

Technical data

Lowest CSP position*	509 mm
Highest CSP position*	1609 mm
Lifting range (range of heights)	1100 mm
Total length	1158 mm
Total width (closed) / (opened)	565 - 1026 mm
Folded length	1236 mm

Folded height	294 mm
Folded width	565 mm
Turning circle	1188 mm
Weight	35 kg
Maximum user weight	130 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.
* CSP = Central suspension point

VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road
HR4 9NS Hereford - United Kingdom
Company registration number 16075613

www.vermeiren.co.uk
info@vermeiren.co.uk



For more information, contact your local dealer

VERMEIREN
we care for you