



Forest KIDS

Powered wheelchair kids

Ready for all terrains! The Forest Kids is an powered wheelchair with an attractive design. The seat can be optimally adjusted for maximum seating comfort for every child.



VERMEIREN
we care for you

All terrains

Some features include: height-adjustable handles, the seat and the back are adjustable thanks to electric springs. All four wheels are always in contact with the ground, which guarantees optimum stability on all terrains.

- Adjustable to suit its users
- Reclinable seat and backrest with powered springs
- Durable upholstery
- Excellent performance in-and outdoor
- Qualitative and trustworthy
- Stability during driving

Colours



Electric-powered reclinable backrest and seat



Front and rear lights for perfect visibility during the darker moments of the day



Anti-tipping and taxi hooks for additional safety

Technical data

Total width	600 mm
Total height	930 - 1080 mm
Total length	1050 mm
Seat height	510 mm
Seat width	320 / 350 / 370 mm
Seat depth	330 - 390 mm

Batteries	2x12V/75Ah/220W/350W
Speed*	6 km/h or 10 km/h
Driving range*	30 km
Weight	95 kg
Maximum user weight	50 kg
Additional info	

* Theoretical range, calculated in accordance with EU standards under ideal conditions. Speeds and range may deviate (approx. 15%). Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road
HR4 9NS Hereford - United Kingdom
Company registration number 16075613

www.vermeiren.co.uk
info@vermeiren.co.uk



Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Forest Kids** complies with the highest norms, laid down by certified test institutes in most European countries.

For more information, contact your local dealer