

# Springer

## Powered Wheelchair Kids

A wheelchair designed to suit the individual needs of each kid. With its many adjustable features, Springer allows for optimal adjustments to ensure ultimate comfort and driving fun.



**VERMEIREN**  
we care for you

## Social integration

Springer is a perfect choice for kids. Easily adapt the chair to any social situations thanks to the option heighten your eye level or to reach high places.

- 300 mm seat lift included
- Electric-powered reclining seat
- LED front and rear lights
- Anti-tipping wheels
- Ground clearance of 80 mm
- Sturdy front and rear wheels with studs on the rear wheels for extra stability



300mm seat lift included

## Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

## Quality guaranteed

Regarding quality, the **Springer** complies with the highest norms, laid down by certified test institutes in most European countries.



## Colours



Electric-powered inclinable seat and backrest

## Technical data

Total width	600 mm
Total height	960 - 1230 mm
Total length	1010 mm
Seat height	530 - 810 mm
Seat width	260 - 380 mm
Seat depth	310 - 370 mm

Batteries	2x12V/50Ah/220W or 2x12V/50Ah/350W
Speed*	6 km/h or 10 km/h
Driving range*	23 km
Weight	90 kg
Maximum user weight	50 kg
Additional info	 

\* Theoretical range, calculated in accordance with EU standards under ideal conditions. \* Speeds and range may deviate (approx. 15%). Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

## VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road  
HR4 9NS Hereford - United Kingdom  
Company registration number 16075613

[www.vermeiren.co.uk](http://www.vermeiren.co.uk)  
[info@vermeiren.co.uk](mailto:info@vermeiren.co.uk)



For more information, contact your local dealer

**VERMEIREN**  
we care for you