

Eclips XXL

XXL wheelchair

Aluminum foldable wheelchair with double cross.
The Eclips XXL provides a natural comfort with its anatomical seat and backrest. This wheelchair, equipped with reinforced rear wheels, can support a maximum weight of 250 kg.

Easy transfers thanks to the foldable aluminum footplate in one piece.



VERMEIREN
we care for you

Key features



Anatomical seat and backrest



Reinforced side plates



Foldable aluminium footplate



Adjustable heel strap

Options



**B83
Canister holder**



**B20/B58
Safety belts**



**B42
Shopping net**



B74 Drum brake



B78 Anti-tipping wheels



B77 Tipping piece

Colours



C27 Blue

Technical data

Total width	850-900-950-1000 mm
Total height	1070 mm
Total length	1140 mm
Folded width	375 mm
Seat width	600-650-700-750 mm
Height armrests	200 mm

Seat height	530 mm
Height backrest	620 mm
Seat depth	530-540-550-560 mm
Weight	38,5 - 40,5 kg
Maximum user weight	250 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

VERMEIREN UK LTD

Suite 1&2 Marshall Business Centre - Faraday Road
HR4 9NS Hereford - United Kingdom
Company registration number 16075613

www.vermeiren.co.uk
info@vermeiren.co.uk



For more information, contact your local dealer

VERMEIREN
we care for you