

VERMEIREN, manufacturer of wheelchairs since 1957, is a major player in the world of medical equipment. Over the years, the Vermeiren Group has specialized in producing wheelchairs, nursing beds, shell chairs, comfort chairs and any other product that may accompany a patient or user in maintaining his autonomy and mobility.



Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.



VERMEIREN UK LTD
Suite 1&2 Marshall Business Centre - Faraday Road - HR4 9NS Hereford
United Kingdom - Company registration number 16075613

VERMEIREN
we care for you

R.E.: VERMEIREN group, Vermeirenplein 1 - 2820 Kalmthout, Belgium - SB Rollators UK 2026-08 v1.0

Altima

PREMIUM ROLLATOR

MOVE LIGHT, LIVE COMFORTABLY

Regain your confidence and enjoy greater independence with the **Altima Rollator**. Designed for comfort, practicality, and ease of use, the Altima is a lightweight, foldable 4-wheel rollator made from durable aluminum, making it easy to maneuver, transport, and store.

For optimum comfort, the Altima is fitted with an extra-comfortable seat and a removable backrest made from soft material. This allows the user to relax for a moment whenever they wish. The elegant, removable shopping bag also offers plenty of space.



british green



matt beige



The solid tyres ensure a pleasant driving experience on a variety of road surfaces, whilst the refined finish and stylish colours - British Green and Matt Beige - further emphasise the Altima's premium look. Even with its convenient closing system, the Altima maintains an

elegant appearance when stored away. The Altima is also fitted with multiple reflectors at the front, rear and sides. Designed for convenience and ease of use, the Altima lets you enjoy safe, comfortable mobility straight away.

Features



removable belt support



removable shopping bag



height-adjustable push handles



safety lock during folding

Technical information

■ Total width	610 mm
■ Total length	700 mm
■ Total height	810 - 985 mm
■ Seat width	425 mm
■ Seat height	640 mm
■ Folded width	245 mm

■ Seat depth	195 mm
■ Turning radius	870 mm
■ Wheels	210 x 36 mm
■ Weight	8,25 kg - 7,35 without basket and backrest
■ Maximum user weight	130 kg
■ Additional information	CE

ROLLATORS



VERMEIREN
we care for you

Four-Light

The Four-light is an elegant and trendy 4-wheel rollator. Available in a sleek carbon grey design and complete with parking brakes, back belt support, a shopping net, and a crutch holder. The handles are adjustable in height and the Four-Light is easily foldable for smooth and convenient transport.



Total width	660 mm
Total length	750 mm
Total height	870 - 990 mm
Seat width	450 mm
Seat height	550 mm
Weight	8,45 kg
Maximum user weight	130 kg

Quava

Quava is a modern, lightweight, foldable rollator with multifunctional brakes. Ideal for the car, travel , errands or for shopping. Thanks to the wide padded back belt you can easily sit down when the brakes are locked.



Total width	655 mm
Total length	710 mm
Total height	830 - 955 mm
Seat width	450 mm
Seat height	540 mm
Weight	7,05 kg
Maximum user weight	136 kg

ECO-Light II

The Eco-Light II is a foldable rollator with 4 wheels, equipped with a removable basket, a padded seat and a crutch holder. A collapsible bar offers additional support while seated. Complete with parking brakes and height-adjustable handles.



Total width	615 mm
Total length	695 mm
Total height	800 - 930 mm
Seat width	355 mm
Seat height	560 mm
Weight	7,5 kg
Maximum user weight	120 kg

286i

Foldable 4 wheel rollator. Equipped with a removable transport basket, a plexi tray-tablet, a seat and a cane holder. Parking brakes and height-adjustable handles. Folding model thanks to a very simple and reliable system.



Total width	600 mm
Total length	680 mm
Total height	780 - 990 mm
Seat width	420 mm
Seat height	560 mm
Weight	9,85 kg
Maximum user weight	130 kg

Quadri-Light

A 4-wheel, folding, lightweight aluminum rollator with brakes is ideal for people with an irregular stride but able to brake when they feel slightly off balance. A seat and a support strap allow the user to rest after long use.



Total width	650 mm
Total length	700 mm
Total height	805 - 930 mm
Seat width	450 mm
Seat height	610 mm
Weight	7,25 kg
Maximum user weight	130 kg

Goliat

The Goliat is a new 4-wheel rollator solid with a maximum user weight of 200 Kg. Equipped with a retaining strap, brakes with parking position and height-adjustable handles. It is possible to adjust its height at the level of the wheels at the front and at the rear.



Total width	650 mm
Total length	750 mm
Total height	830 - 1020 mm
Seat height	490 - 560 mm
Seat depth	410 mm
Weight	10,7 kg
Maximum user weight	200 kg

ECO-Plus

The Rollator ECO-Plus is a 4-wheel rollator with the feature of having an adjustment on the front and the back allowing to adjust the height of his seat. Equipped with a net, a foam cushion, cane holder and a support belt provides additional comfort.



Total width	610 mm
Total length	680 - 795 mm
Total height	845 - 1020 mm
Seat width	360 mm
Seat height	555 - 625 mm
Weight	7,5 kg
Maximum user weight	136 kg

Goliat Artro

Heavy-duty upright rollator with forearm supports. Goliat Artro has a modern and strong design with height-adjustable forearm supports to relieve stress on hands and wrists whilst standing upright. A maximum user weight of 200 kg makes this rollator perfect for users who need extra stability.



Total width	650 mm
Total length	750 mm
Total height	1140 - 1310 mm
Seat height	490 - 560 mm
Seat depth	410 mm
Weight	12,35 kg
Maximum user weight	200 kg