




Bill of Material			
No.	Plan	Article	Quantity
1	B31 Plastic	1906616	1
2	P00083426	2691154	1
3	Nut M6 DIN985	2994063	2
4	Washer M6 DIN125	2994414	3
5	M6x35 DIN7984	2994654	2
6	Cap D25	1991857	1
7	M6x8 ISO7380	2994901	1
8	XP02043	2993671	1
9	M4x12 DIN7984	1905782	1
10	M4 DIN985	2994901	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
 EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
 FROM VERMEIREN GROUP.

1921190

B31 CRUTCH CARRIER V500 LIGHT

Plan No.

E00131562/B1.1

Scale: 1:4	Drafted	Verified	Confirmed
Date	07/07/2026	24/07/2026	24/07/2026
Initials	PM	PM	AZ

