








Detail A
Scale: 1:5



Detail B
Scale: 1:5

④ Note: Serial number sticker must be fixed parallel to frame.



Detail C
Scale: 1:5

Material				
No.	Title	Article	Sticker	Quantity
1	Sticker Ø logo Vermeiren	T001002		2
2	Sticker Vermeiren 85x15	T001051		2
3	Sticker taxi hook	T001076		4
4	Sticker serial number	T001097		1

STICKER CONFIGURATION JAZZS50

STICKER CONFIGURATIE JAZZS50

Plan No.

E001127/C

Scale: 3:40

Drafted

Verified

Confirmed

Date

01/04/2019

16/07/2026

16/07/2026

Initials

TT

MP

MP



E001127 / C