


BK BLACK
RD RED
LBL LIGHT BLUE
BL BLUE
LGR LIGHT GREEN
GR GREEN
VI VILLET
OR ORANGE
YE YELLOW
WH WHITE
LRD LIGHT RED
PK PINK
GY GRAY
/ AND
I OR



#	Reference	Art number	Description	Qty	OPTION
1	EC_BAT1, EC_BAT2	1917919	EC BAT 80Ah Carbon	2	80Ah
		1916507	EC BAT 80Ah Gel	2	PL 80Ah
		1917300	EC BAT 100Ah AGM	2	100Ah
2	EC_COM1	1923231	E BATTERIE 12V DC EVX12-100 GEL LANDPORT	2	PL 100Ah
		1922287	SV KUARZ KEY SETS BH 2801	1	KEY01
		1922288	SV KUARZ KEY SETS BH 2802	1	KEY02
		1922289	SV KUARZ KEY SETS BH 2803	1	KEY03
		1922290	SV KUARZ KEY SETS BH 2804	1	KEY04
		1922291	SV KUARZ KEY SETS BH 2805	1	KEY05
		1922292	SV KUARZ KEY SETS BH 2806	1	KEY06
		1922293	SV KUARZ KEY SETS BH 2807	1	KEY07
		1922294	SV KUARZ KEY SETS BH 2808	1	KEY08
		1922295	SV KUARZ KEY SETS BH 2809	1	KEY09
3	EC_COM2	3992019	EC COM TACOS 5K VR + ALU BLOCK BH I	1	
4	EC_COM3	3992014	EC COM SV CARPO2 BUZZER BH	1	
5	EC_COM4	3992016	EC COM SV CARPO2 USB BH	1	
6	EC_COM5	1921871	EC COM SV KUARZ OPTION CAMERA FOR PCB	1	CAM
7	EC_COM6	1921864	EC COM KUARZ PCB	1	
8	EC_COM7	1921861	EC COM SV KUARZ SPEED SENSOR	1	
9	EC_COM8	1921863	EC COM SV KUARZ CHARGE HOLDER	1	
10	EC_COM9	3990624	EC COM KUARZ CONTROLLER S DRIVE 140	1	SE XD S
		3990623	EC COM KUARZ CONTROLLER S DRIVE 120	1	SE
11	EC_L2	1920893	EC L SV KUARZ FRONT LIGHT	1	
12	EC_L1, EC_L3	1920894	EC L SV KUARZ/SYNC/EON FRONT BLINKER	2	
13	EC_L4, EC_L7	1920895	EC L SV KUARZ/SYNC/EON REAR BLINKER	2	
14	EC_L5, EC_L6, EC_L8, EC_L9	1920896	EC L SV KUARZ/SYNC/EON REAR/BRAKE LIGHT	4	
15	EC_MOT1	3993315	EC MOT KUARZ (SE) MOTOR 750W 1:20	1	SE
		3991650	EC MOT KUARZ SV CARPO2 MOTOR 750W XD COMPLETE SET 1:22	1	XD
16	EC_SW1	3992376	EC SW BATTERY SWITCH 150A A701	1	
18	EC_WIRE2	1921858	EC WIRE KUARZ CABLE FOR FRONT LIGHTS	1	
20	EC_WIRE4	1922286	EC WIRE SV KUARZ MAIN FUSE 175A	1	
21	EC_WIRE5	1921860	EC WIRE SV KUARZ CONTROL CABLE	1	
22	EC_WIRE6	1921862	EC WIRE KUARZ MAIN POWER CABLE	1	
23	EC_WIRE7, EC_WIRE8	1921859	EC WIRE SV KUARZ CABLE FOR REAR LIGHTS	2	
24	UNDEPICTED	2998581	CABLE TIES AHM PM-150	13	
25	UNDEPICTED	1993111	BINDING BAND NYLON L=500MM	6	
26	UNDEPICTED	3994914	FERRITE CORE KCF-130-B PK2503CM	1	
27	UNDEPICTED	3994915	FERRITE CORE K5B-RC-16X28X9-M2 PK1602IM	2	CAM
28	UNDEPICTED	3994916	FERRITE CORE K5B-RC-20X30X13-MBP2007IMS	1	CAM
29	UNDEPICTED	1923012	BATTERY IVRA C-GO SMART 8A AGM/GEL	1	80Ah
30	UNDEPICTED	1923011	BATTERY IVRA C-GO SMART 10A AGM	1	100Ah
31	UNDEPICTED	1916118	CHARGER VICTRON IMPULSE SMART 8A AGM/GEL	1	CH 80Ah; PL 80Ah; PL 100Ah

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORIZATION FROM VERMEIREN GROUP

BOM EC KUARZ

Size: A4	Drafted	Verified	Confirmed
Date:	2026-06-24	2026-06-24	2026-06-24
Initials:	FM	PM	MK

2998640/C

Plan No.

2998640/C

