




NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
 EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
 FROM VERMEIREN GROUP.

FOREST 3+ FRONT LOWER FRAME			
Scale:	1:4	Drafted	Verified
Date	22.11.2016	23/06/2026	23/06/2026
Initials	SD	KA	AZ

2993393

Plan No.

E001057/G



Material			
No.	Plan	Article	Quantity
1	M8	1902024	2
2	X001119	1001554	8
3	X001772	1521500	2
4	X005141	1910010	1
5	X005589	1912549	2
6	X005614	1912738	2
7	Y018567	17734..	2
8	Y022165	20367..	1
9	---	1001409	2
10	M12 type M	1911902	1
11	M12 DIN125	1001561	2
12	M6 DIN125	1001564	2
13	D38	1001611	2
14	M12 DIN985	1001684	3
15	M6x50 DIN912	1007970	2
16	/	/	/
17	M6	1008492	2
18	M12 D16x80 ISO7379	1907488	1
19	M8x120 DIN912	1994676	2
20	Black Line R-Fill	1915162	1
21	Grey Block Air	1915158	1
	Black Line Air	1915160	1
	Black Line R-Fill	1915163	1
	Grey Block Air	1915158	1
	Black Line Air	1915161	1
22	M8x12 DIN912	1008963	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
 EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
 FROM VERMEIREN GROUP.

2993393

FOREST 3+ FRONT LOWER FRAME				Plan No.
Scale: 1:4				E001057/G
Date	22.11.2016	23/06/2026	23/06/2026	
Initials	SD	KA	AZ	