




No.	Plan	Article	Qty
1	E00127220	-	1
2	P00118924 DP3	2920054	1
3	P00118924 DP4	2920954	1
4	P00116672	GP331	1
5	P00116753	GP335	1
6	P00118936	GP336	2
7	DIN7984 M8X75 BLACK	2994912	3
8	1915949 M8x90	1915949	1
9	P00118829	GP332	1
10	WASHER M8 D17XD9X1.6 DIN126 BLACK ZINC	1996214	8
11	1001556WasherM8DIN125B	1001556	4
12	WASHER PVC M12 DIN433 D13X20X2	2993281	6
13	SELFLOCKING NUT M8 DIN985 BLACK	1996213	4
14	2994067 M6x16	2994067	2
15	PLASTIC WHEEL D55XD8	1906057	2
16	1001517WingedNutM8	1001517	1
17	M8x30 DIN933/8.8	1001443	1
18	2994415 WASHER PA6 M8 DIN125	2994415	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

LIFT MECHANISM BLUMIL SPORT

Scale: 1:5	Drafted	Verified	Confirmed
Date	27.05.2026	27.05.2026	27.05.2026
Initials	KK	PM	MK

Plan No.

E00118924/A.1

