





NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 1/7

PLEGO

Plan No.


E00052670/I.1

Scale: 1:10	Drafted	Verified	Confirmed
Date	27/09/2023	09/06/2026	09/06/2026
Initials	BK	DR	AS



* Use WEICONLOCK AN 302-43

** Optional equipment





* Use:
WEICONLOCK AN 302-43

Sheet 2/7

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

PLEGO				---
				Plan No.
				E00052670/I.1
Scale:	1:10	Drafted	Verified	Confirmed
Date	27/09/2023	09/06/2026	09/06/2026	
Initials	BK	DR	AS	







NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
 EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
 FROM VERMEIREN GROUP.

Sheet 3/7

PLEGO				---
				Plan No.
				E00052670/I.1
Scale:	1:10	Drafted	Verified	Confirmed
Date	27/09/2023	09/06/2026	09/06/2026	
Initials	BK	DR	AS	







* Use:
WEICONLOCK AN 302-43

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

Sheet 4/7

PLEGO				---
				Plan No.
				E00052670/I.1
Scale: 1:10	Drafted	Verified	Confirmed	
Date	27/09/2023	09/06/2026	09/06/2026	
Initials	BK	DR	AS	




* Use:
WEICONLOCK AN 302-43


NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

Sheet 5/7

PLEGO				---
Scale: 1:10				Plan No.
Date				E00052670/I.1
Initials				
Drafted	27/09/2023	Verified	09/06/2026	
Confirmed	09/06/2026	AS		



FOR STICKERS T001208-T001216
USE GAUGE P00068642



NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

Sheet 6/7

PLEGO			
Scale: 1:10	Drafted	Verified	Confirmed
Date	27/09/2023	09/06/2026	09/06/2026
Initials	BK	DR	AS

Plan No.
E00052670/I.1

No.	Plan	Article	Qty	24	WASHER M6	2994414	16	45	SCREW PARKER 4.8X19 DIN7983	1001534	10								
1	E00063648	---	1	25	WASHER M10 DIN125	2998106	4					46	A00056460	2997888**	1				
2	E00062349	---	1*					26	M8x25	1990501	2					47	A00056459	2997885	1
2**	E00119881	---	1*	27	NUT M8	2995398	2	48	M5x25	2994083	2								
3	E00062350	---	1	28	M6X25 DIN912	2996887	8	49	NUT M5	2994006	2								
4	E00062352	---	1	29	X000862	1006463	8	50	M6X10	1904683	4								
5	E00056068	---	1	30	NUT DIN985 M6	2994063	4	51	E00054750	---**	1								
6	E00063644	---	1	31	M6x45	2994086	4	52	E00065217/ E00066678	---**/ ---**	1								
7	E00062355	---	1	32	P00053220	26742B8 26742A9 26742B2 26742A7 26742D5 2674221	8					53	P00066352	T001061	2				
8	P00053716	26771	2					54	P00066353	T001060	2								
9	P00062160	2677665	1					55	A00066354	1920500**	1								
10	P00061760	2677765	1					56	E00066356	---**	1								
11	P00053715	2678965	1					57	E00066564	---**	1								
12	Battery 24V 12Ah Battery 24V 17Ah	1918834 1919283**	1 1					58	E00066475	---**	1								
13	XC01046	1919584	2	33	M4X10 DIN7991	2998006	8	59	A00066482	2998031**	1								
14	M6x16	2994067	4	34	X000106	G003218	4	60	P00066485	1920046**	2								
15	NUT M6	2994063	4	35	M8 DIN125	2994415	8	61	NUT M8	1996213	2								
16	-	-	-	36	M8X60 ISO7380	2995491	2	62**	B31 CRUTCH HOLDER 25mm	1906616	1								
17	1919240	1919240	1	37	M5x12	2994074	4					63	Drive train standard Drive train + SE19 lights**	2998597 2998582	1 1				
18	2993352 ELEMENT 1 /2993352 ELEMENT 2	2993352	1	38	P00054831	26779	2	64**	Sticker "SPRING"	T001285	2								
	39			A00056456	2997887	1	65					Sticker Serial Number	T001097	1					
	19			P00061690	T001215 T001214 T001212 T001216 T001209 T001208	2		40	A00056458	2997886	1								
41		A00056457	2997884				1												
42		E00052668	---				1												
43		STAR KNOB M6x35 VGF	2990215				1												
20	P00053212	2997848	2	44	X001473	1002994	10												
21	M8 DIN125B	1001556	2																
22	M6x20	2994004	4																
23	P00053209	2997847	4																
												* 1 pair contains 2 pcs - check AS system ** Optional equipment				NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.			
								PLEGO				Plan No.							
												E00052670/I.1							
								Scale: 1:10				Drafted		Verified		Confirmed			
								Date				27/09/2023		09/06/2026		09/06/2026			
								Initials				BK		DR		AS			