



No.	Plan	Article	Qty
1	P00119285	2921954	2
2	P00116717	2920354	2
3	P00116709	2920254	2
4	P00118936	GP336	1
5	P00116672	GP331	1
6	P00115681	GP333	2
7	P00116714	GP334	2
8	P00120980	2921154	1
9	P00121226	2980011	1
10	2994127BoltM6x10black	2994127	1
11	1004414SpringLockWasherDIN127M6	1004414	1
12	1915949 M8x90	1915949	1
13	DIN7984 M8X75 BLACK	2994912	1
14	1005289InnerHexM8x70	1005289	1
15	M8 D17XD9X1.6 DIN126 BLACK ZINC	1996214	6
16	WASHER PVC M12 DIN43	2993281	4
17	SELFLOCKING NUT M8 DIN985 BLACK	1996213	3
18	M6X10 DIN912/8.8	1001491	2
19	DIN125A M6 PA6	1000973	2
20	MOTOR BLUM.GO 36VDC 80W	2980063	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

## LIFT MECHANISM 1ST STEP BLUMIL GO

Plan No.  
E00127223/A.1

Scale:	1:5	Drafted	Verified	Confirmed
Date	27.05.2026	27.05.2026	27.05.2026	27.05.2026
Initials	KK	PM	MK	

