





Detail A  
Scale: 1:3

Drawing shows only right side

Notes:

\*use only with 2993810 shortened

EP00266 / J

| BACKREST TRIGO |            |            |            | Plan No.                                                                              |
|----------------|------------|------------|------------|---------------------------------------------------------------------------------------|
| Page1/3        |            |            |            | EP00266/J                                                                             |
| Scale: 1:7     | Drafted    | Verified   | Confirmed  |  |
| Date           | 06/11/2017 | 27/07/2026 | 27/07/2026 |                                                                                       |
| Initials       | AS         | DZ         | AS         |                                                                                       |

| Material (25-30) |         |         |          |
|------------------|---------|---------|----------|
| No.              | Plan    | Article | Quantity |
| 3                | YP05930 | 2785565 | 1        |
| 4                | YP06265 | 2785265 | 1        |

| Material (30-35) |           |         |          |
|------------------|-----------|---------|----------|
| No.              | Plan      | Article | Quantity |
| 3                | YP05930   | 2785565 | 1        |
| 4                | P00025189 | 2643065 | 1        |

| Material (35-40) |           |         |          |
|------------------|-----------|---------|----------|
| No.              | Plan      | Article | Quantity |
| 3                | YP05932   | 2786065 | 1        |
| 4                | P00025189 | 2643065 | 1        |

| Material (40-45) |           |         |          |
|------------------|-----------|---------|----------|
| No.              | Plan      | Article | Quantity |
| 3                | YP05933   | 2786165 | 1        |
| 4                | P00025189 | 2643065 | 1        |

| Material (45-50) |         |         |          |
|------------------|---------|---------|----------|
| No.              | Plan    | Article | Quantity |
| 3                | YP05932 | 2786065 | 1        |
| 4                | YP06269 | 2785465 | 1        |

| Material |              |         |          |
|----------|--------------|---------|----------|
| No.      | Plan         | Article | Quantity |
| 5        | YP05910 (25) | 2786465 | 1        |
|          | YP05911 (27) | 2786565 | 1        |
|          | YP05912 (29) | 2786665 | 1        |
|          | YP05913 (31) | 2786765 | 1        |
|          | YP05914 (33) | 2786865 | 1        |
|          | YP05915 (35) | 2786965 | 1        |
|          | YP05916 (37) | 2787065 | 1        |
|          | YP05917 (39) | 2787165 | 1        |
|          | YP05918 (41) | 2787265 | 1        |
|          | YP05919 (43) | 2787365 | 1        |
|          | YP05920 (45) | 2787465 | 1        |
|          | YP05921 (47) | 2787565 | 1        |
|          | YP05922 (49) | 2787665 | 1        |
|          | YP05906 (51) | 2787765 | 1        |
| 6        | YP05995 (25) | 2789465 | 1        |
|          | YP05997 (27) | 2789565 | 1        |
|          | YP05998 (29) | 2789665 | 1        |
|          | YP05999 (31) | 2789765 | 1        |
|          | YP06000 (33) | 2789865 | 1        |
|          | YP06001 (35) | 2789965 | 1        |
|          | YP05992 (37) | 2790165 | 1        |
|          | YP06002 (39) | 2790065 | 1        |
|          | YP06003 (41) | 2790265 | 1        |
|          | YP06004 (43) | 2790365 | 1        |

| Material |                   |         |          |
|----------|-------------------|---------|----------|
| No.      | Plan              | Article | Quantity |
| 6        | YP06005 (45)      | 2790465 | 1        |
|          | YP06006 (47)      | 2790565 | 1        |
|          | YP06007 (49)      | 2790665 | 1        |
|          | YP06008 (51)      | 2790765 | 1        |
| 7        | YP06896 (25)      | 2765765 | 1        |
|          | YP06897 (27)      | 2765865 | 1        |
|          | YP06898 (29)      | 2765965 | 1        |
|          | YP06899 (31)      | 2766065 | 1        |
|          | YP06900 (33)      | 2766165 | 1        |
|          | YP06901 (35)      | 2766265 | 1        |
|          | YP06902 (37)      | 2766365 | 1        |
|          | YP06903 (39)      | 2766465 | 1        |
|          | YP06904 (41)      | 2766565 | 1        |
|          | YP06905 (43)      | 2766665 | 1        |
| 8        | YP06906 (45)      | 2766765 | 1        |
|          | YP06907 (47)      | 2766865 | 1        |
|          | YP06908 (49)      | 2766965 | 1        |
|          | YP06909 (51)      | 2767065 | 1        |
| 8        | Fold. push handle | EP00246 | 1        |
|          | Y000462           | F000050 | 1        |
| 9        | Y004353           | F000411 | 1        |
| 10       | XP01560           | 2993561 | 4        |
| 11       | X000273           | G0606   | 1        |

\*\* When using a folded handle, select 2 sizes shorter tube.

EP00266 / J

| BACKREST TRIGO |            |            |            | Plan No.                                                                              |
|----------------|------------|------------|------------|---------------------------------------------------------------------------------------|
| Page 2/3       |            |            |            | EP00266/ J                                                                            |
| Scale: 1:7     | Drafted    | Verified   | Confirmed  |  |
| Date           | 06/11/2017 | 27/07/2026 | 27/07/2026 |                                                                                       |
| Initials       | AS         | DZ         | AS         |                                                                                       |

| Material |                 |         |          |
|----------|-----------------|---------|----------|
| No.      | Plan            | Article | Quantity |
| 13       | M6              | 2994414 | 2        |
| 14       | M6x45           | 2994086 | 2        |
|          | M6x55*          | 2994376 | 2        |
| 15       | ---             | ---     | ---      |
| 16       | ---             | ---     | ---      |
| 17       | CAP D16         | 1001598 | 1        |
| 18       | M6x10           | 2994127 | 4        |
| 19       | Hand grip GR-KS | 2993809 | 1        |
|          | GR-ERGO-S       | 2993810 | 1        |
| 20       | M6              | 2994063 | 1        |
| 21       | YP06168(39)     | 2797065 | 1        |
|          | YP07800(41)     | 2782265 | 1        |
|          | YP07801(43)     | 2782365 | 1        |
|          | YP07802(45)     | 2782465 | 1        |
|          | YP07803(47)     | 2782565 | 1        |
|          | YP07804(49)     | 2782665 | 1        |
|          | YP07863(51)     | 2600665 | 1        |
| 22       | Y017346         | 1748054 | 1        |
| 23       | Y017412         | 1747800 | 1        |
| 24       | Ø3x16 DIN1481   | 1003791 | 1        |
| 25       | M6x15           | 1906622 | 1        |
| 26       | X004876         | 17477   | 1        |
| 27       | Handgrip        | 1906244 | 1        |

\* Use with B82 fixed, B82S fixed and B82S C fixed

EP00266 / J

| BACKREST TRIGO |            |            |            | Plan No.                                                                              |
|----------------|------------|------------|------------|---------------------------------------------------------------------------------------|
| Page 3/3       |            |            |            | EP00266/ J                                                                            |
| Scale: 1:7     | Drafted    | Verified   | Confirmed  |  |
| Date           | 06/11/2017 | 27/07/2026 | 27/07/2026 |                                                                                       |
| Initials       | AS         | DZ         | AS         |                                                                                       |