




No.	Plan	Article	Qty
1	BRACKET LEFT	2928804	1
2	BRACKET RIGHT	2928904	1
3	NUT M8 DIN985	3991743	2
4	WASHER M8	1003199	2
5	SCREW M8X20	1992349	2
6	SEAT BELT	2993352	1



NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
FROM VERMEIREN GROUP.

SV KUARZ SEAT BELT COMPL.				---
				Plan No. E00130927/A.1
Scale: 1:10	Drafted	Verified	Confirmed	
Date	29/06/2026	29/06/2026	29/06/2026	
Initials	BK	BK	MK	