



No.	Plan	Article	Qty
1	E00075734	---	1
2	P00049512	3991104	1
3	P00049511	3990780	1
4	T001090	T001090	1
5*	P00019685	T001069/T001068	1
6	T001263	T001263	1
7	T001261	T001261	1
9	T001280	T001280	1
10	T000054	T000054	1


* Only for the Swiss market
 JIG P0122198 for position 6, 7, 8, 9
 JIG P0121969 for certification label
 T000054 - Sticker on box

Certification label area


Stick on inner wall of Battery Cover
 in showed place



A (1:5)



B (1:7)



NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
 EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
 FROM VERMEIREN GROUP.

SV KUARZ STICKER COFIGURATION				Plan No. E00122492/C.1
Scale: 1:15	Drafted	Verified	Confirmed	
Date	19/02/2026	30/06/2026	30/06/2026	
Initials	BK	BK	MK	