




No.	Plan	Article	Qty
1	P00114423	2919404	1
2	TUBE PLUG 30x30 Black	1003198	6
3	SCREW M8X45 DIN912	2995397	4
4	Star Knob VGF 40 M8x20	1001580	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
 EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
 FROM VERMEIREN GROUP.

SV KUARZ OPTIONAL EQUIPMENT BAR			
Scale: 1:3	Drafted	Verified	Confirmed
Date	13/11/2025	23/03/2026	23/03/2026
Initials	BK	BK	MK

Plan No.
E00116546/B.1