










Front View



Right View



B (3:20)



Note: Serial number sticker must be fixed parallel to frame.

Material			
No.	Title	Article	Quantity
1	Sticker Dia 35 logo Vermeiren	T001002	2
2	Sticker serial number	T001097	1
3	Sticker Bobby EVO white	T001252	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

STICKER CONFIGURATION BOBBY EVO

Scale: 3:40	Drafted	Verified	Confirmed
Date	23-06-2026	23-06-2026	23-06-2026
Initials	DT	MP	MP

Plan No.

A00130189/A

