




Material							
No.	Plan	Article	Quantity	No.	Plan	Article	Quantity
1	X001016	1001682	4	18	M10 DIN985	1001683	2
2	X002282	1992854	4	19	M8x45 DIN912	1003110	4
3	X002792	1555954	4	20	M6x12 ISO7380	1007000	9
4	X003403	1617400	4	21	M5x35 DIN912	1007961	2
5	X004778	1908255	4	22	M12 DIN125	1993362	2
6	X005576	1912735	1	23	---	1906791	2
7	X005609	1912734	1	24	---	1906792	2
8	X005637	1912761	1	25	---	1908165	2
9	X005638	1912762	1	26	M5x10	1909204	2
10	Y020969 DP	20363..	1	27	80A	1909646	2
11	Y021565	2025854	1	28	150A	1910488	1
12	Y022051	2025954	1	29*	---	1912309	1
13	M6 DIN125	1000973	6	30	LINX 120A	1912415	1
15	M5 DIN125	1001550	2	31	M6	1912736	3
16	D28	1001606	2	32	M10x65 DIN912	1912777	2
17	M5 DIN985	1001679	2	33	---	1999843	2
				34	M5x10 ISO7380	2991862	4
				35	P00126766	2980107	1

* Option SE50

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
 EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION
 FROM VERMEIREN GROUP.

Sheet 2/2

2993391

Forest 3+ Rear Lower Frame

Plan No.

E001058/E

Scale: 1:5	Drafted	Verified	Confirmed
Date	23.11.2016	30/07/2026	30/07/2026
Initials	SD	AZ	MK

