


NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 1/2

| | | | | |
|--|------------|------------|------------|---|
| Electric backrest frame Springer Kids | | | | Plan No. EP00029/G |
| Scale: 1:1 | Drafted | Verified | Confirmed |  |
| Date | 23/03/2011 | 11/03/2026 | 11/03/2026 | |
| Initials | AS | PM | MK | |

| Bill of Material | | | |
|------------------|---------------|----------|----------|
| No. | Plan | Article | Quantity |
| 1 | YP02544 | 2803354 | 1 |
| 2 | YP02540 | 2803254 | 1 |
| 3 | YP03008 | 28222... | 1 |
| 4 | --- | 1003299 | 1,2m |
| 5 | --- | 1906790 | 1 |
| 6 | YP02160 | 1893554 | 2 |
| 7 | YP02159 | 1893454 | 2 |
| 8 | XP00550 | 2990635 | 2 |
| 9 | M6x12 | 1007000 | 4 |
| 10 | M6X16 | 1008459 | 4 |
| 11 | DIN7991-M6x35 | 2991297 | 4 |
| 12 | DIN125A-M6 | 1000973 | 4 |
| 13 | DIN985-M6 | 1001680 | 4 |
| 14 | DIN7991-M8x45 | 2991298 | 2 |
| 15 | M10x50 | 1009335 | 1 |
| 16 | DIN125A-M10 | 1001558 | 2 |
| 17 | DIN985-M10 | 1001683 | 1 |
| 18 | X000862 | 1006463 | 2 |
| 19 | DIN985-M8 | 1001682 | 2 |
| 20 | M8x15 | 1006661 | 1 |
| 21 | Blue | 2991367 | 1 |
| | Orange | 2991368 | 1 |
| 22 | M10x40 | 1991481 | 1 |

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 2/2

Electric backrest frame Springer Kids

Plan No.

EP00029/G

| | | | | |
|----------|------------|------------|------------|-----------|
| Scale: | 1:1 | Drafted | Verified | Confirmed |
| Date | 23/03/2011 | 11/03/2026 | 11/03/2026 | |
| Initials | AS | PM | MK | |

